

Vantive Product Security Bulletin

RE: Mirth Connect Vulnerability, CVE-2023-43208

Publication Date: January 31, 2024

ISSUE SUMMARY

We take cybersecurity at Vantive seriously and are currently monitoring and evaluating the recently published announcement of the vulnerability affecting Mirth Connect under CVE-2023-43208. Mirth Connect is an open-source data integration platform from NextGen HealthCare that is used to connect data from medical systems to an electronic medical record. This vulnerability has been classified as critical. For a more detailed description of this vulnerability, it is recommended customers view the information provided in the [National Vulnerabilities Database](#).

Products in Scope

The products listed below were identified as being impacted by CVE-2023-43208:

- Sharesource Connect

RESPONSE

Customers can upgrade Mith Connect to version 4.4.2 on their Sharesource Connect servers and can contact Vantive Technical Support if any assistance is needed.

If you observe any activity associated with this vulnerability, contact your service representative immediately. Please contact the Vantive Product Security team at global.corp.product.security@vantive.com if you have any additional questions. As further information becomes available, Vantive will assess for any potential impact on its products and update this security bulletin.